


Optional accessories


No.	Description
377-091D	Ring-type unit connecting ring (for MSM-412)
377-092	Ring-type unit connecting ring (for MSM-Z414)



Accessories for Stereo Microscope

Ring Light

- The ring light provides additional, intensive, fluorescent illumination. This type of light source prevents colour errors in the image caused by the colour characteristics of the illumination. The ring-shaped positioning of the light source around the objective ensures even illumination of the subject without interfering shadows.



377-090-1DD

No.	Description
377-090-1DD	Ring-type fluorescent illumination unit
377-090-1E*	Ring-type fluorescent illumination unit (UK only)


Universal mounting support

- For greater convenience when observing large objects that cannot be installed on a standard support.
- The adjustable arm provides for several observation angles.
- Maximum height of object: approx. 225/205mm (MSM-Z414, MSM-465).
- Optimum weight distribution guarantees excellent stability of the instrument.

No.	Applicable Microscope
377-007D*	MSM-412 and MSM-Z414
377-107D*	MSM-465

CCD Camera Adaptor

- 1/2" CCD camera support (camera with load transfer device).
- Switching between binocular and video observation modes.
- Digital diaphragm opening/closing control.
- Under certain observation conditions, visible artifacts may appear in the field of view.



No.	Applicable Microscope
377-094D	MSM-465